

# IMAGEREADY™

## MR-Conditional Pacing System

### Cardiology Order Form

Patient name: \_\_\_\_\_ Date of Birth: \_\_\_\_\_

Pacemaker Model: \_\_\_\_\_ RV Lead Model: \_\_\_\_\_ Atrial Lead Model: \_\_\_\_\_

**According to Boston Scientific's device labeling, the following Conditions of Use must be met in order for your patient to undergo an MRI scan:**

1. Please confirm that your patient has a complete ImageReady MR-Conditional Pacing System, defined as a Boston Scientific MR-Conditional pulse generator and INGEVITY™ MRI Pacing Lead, with all ports occupied by a lead or port plug.

**Resources:** • Boston Scientific MRI Technical Guide, ImageReady MR-Conditional Pacing System • [www.BostonScientific.com/imageready](http://www.BostonScientific.com/imageready)

☐ YES, my patient has a complete ImageReady MR-Conditional Pacing System in the left or right pectoral region and at least six (6) weeks have elapsed since the implantation and/or lead revision or surgical modification of the MR-Conditional Pacing System.

☐ NO, my patient does not have a complete ImageReady MR-Conditional Pacing System.

2. Does your patient have any cardiac-related implanted devices, components, or accessories present other than an ImageReady MR-Conditional Pacing System?

☐ Yes ☐ No

3. Does your patient have any abandoned leads or pulse generators?

☐ Yes ☐ No

4. Does your patient have evidence of a fractured lead or compromised pulse generator-lead system integrity (lead impedances are not within normal range)?

☐ Yes ☐ No

5. Does your patient have an elevated body temperature or compromised thermoregulation at the time of the scan?

☐ Yes ☐ No

6. If your patient is pace-dependent, is the pacing threshold  $\leq 2.0$  V?

☐ Yes ☐ No ☐ NA

7. Before the scan, your patient's pacemaker will be placed in an MRI Protection Mode. How would you like your patient's pacemaker to be programmed?

☐ DOO (Bipolar) Pacing rate: \_\_\_\_\_ ppm ☐ AOO (Bipolar) Pacing rate: \_\_\_\_\_ ppm

☐ VOO (Bipolar) Pacing rate: \_\_\_\_\_ ppm ☐ Pacing Off

☐ MRI Protection Time-out: 12 hours \_\_\_\_\_, 24 hours (nominal) \_\_\_\_\_, 48 hours \_\_\_\_\_ Off \_\_\_\_\_

*MRI Protection Time-out is designed to make sure patients will return to their original pacemaker settings after the scan, without the need of reprogramming.*

- When the Time-out parameter is programmed to a value other than Off, the patient must be out of the scanner before the time programmed elapses
- Warning: If the MRI Protection Time-out value of Off is combined with a Pacing Mode of Off, the patient will not receive pacing until the pulse generator is programmed using a programmer out of MRI Protection Mode and back to normal operation

Physician Signature: \_\_\_\_\_

Physician name: \_\_\_\_\_

Date: \_\_\_\_\_

Clinic Name: \_\_\_\_\_

Clinic Address: \_\_\_\_\_

Clinic Phone: \_\_\_\_\_

*This form may contain patient confidential information. If you receive this form in error, please do not forward it and contact Boston Scientific Technical Services at (800) 227-3422.*

Registered trademarks are the property of their respective owner.

© 2016 Boston Scientific Corporation or its affiliates. All rights reserved. CRM-330611-AC NOV2016