

Abdominal Aorta Ultrasound Protocol

- Proximal Aorta
 - Sag (Gray Scale, Color Doppler, Pulsed Wave)
 - Measurement outer to outer wall
 - Measure peak velocity
 - Trans (Gray Scale)
 - Measurement outer to outer wall
- Mid Aorta
 - Sag (Gray Scale, Color Doppler, Pulsed Wave)
 - Measurement outer to outer wall
 - Measure peak velocity
 - Trans (Gray Scale)
 - Measurement outer to outer wall
- Distal Aorta
 - Sag (Gray Scale, Color Doppler, Pulsed Wave)
 - Measurement outer to outer wall
 - Measure peak velocity
 - Trans (Gray Scale)
 - Measurement outer to outer wall
- Aorta Bifurcation
 - Trans (Gray Scale and Color Doppler)
 - Measurement outer to outer wall RT and LT Iliac
- Right Iliac Artery
 - Sag (Gray Scale, Color Doppler, Pulsed Wave)
 - Measurement outer to outer wall
 - Measure peak velocity
- Right Iliac Artery
 - Sag (Gray Scale, Color Doppler, Pulsed Wave)
 - Measurement outer to outer wall
 - Measure peak velocity